




## Certificate of Analysis

CP-Pre/Adi-62

AxyGen's ready to use buffers, precipitants, salts, additives, and detergents provide the most complete and economical solution to crystallization of proteins, peptides, small molecules and nucleic acids. They can be used as a rapid Initial screening method or to accelerate the creation of optimization/follow up screens. Reagents are supplied in sterile bottles (200ml/bottle).

<b>Product Name:</b>	<b>4.0M Na-Acetate</b>
<b>Description:</b>	Sodium Acetate dissolved in water
<b>Catalog #:</b>	CP-Pre/Adi-62
<b>Chemical Name:</b>	Acetic acid sodium salt
<b>Molecular Formula:</b>	CH <sub>3</sub> COONa
<b>Molecular Weight:</b>	82.03
<b>CAS Number:</b>	<a href="#">127-09-3</a>
<b>Beilstein Registry:</b>	3595639
<b>EC Number:</b>	2048238
<b>Appearance:</b>	White Powder
<b>Purity:</b>	≥99%
<b>Solubility:</b>	Dissolves at 100mg/ml of water yielding a clear, colorless solution
<b>pKa :</b>	4.76 (acetic acid)
<b>MP:</b>	>300 °C



**AxyGen BioSciences**  
33210 Central Ave., Union City, CA 94587  
Ph.: (510) 476-1090 Fax: (510) 476-1091  
Email: [sales@axygenbio.com](mailto:sales@axygenbio.com)  
Website: [www.axygenbio.com](http://www.axygenbio.com)