




Certificate of Analysis

CP-Buffer-20

AxyGen's ready to use buffers, precipitants, salts, additives, and detergents provide the most complete and economical solution to crystallization of proteins, peptides, small molecules and nucleic acids. They can be used as a rapid Initial screening method or to accelerate the creation of optimization/follow up screens. Reagents are supplied in sterile bottles (100ml/bottle).

Product Name:	1.0 M MES, pH 5.5
Description:	MES titrated with NaOH
Catalog #:	CP-Buffer-20
Chemical Name:	2-(N-Morpholino)ethanesulfonic acid
Molecular Formula:	C ₆ H ₁₃ NO ₄ S
Molecular Weight:	195.24
CAS Number:	4432-31-9
Beilstein Registry:	141319
EC Number:	2246323
Purity:	≥99.0%
Appearance:	White powder
Solubility:	Dissolves at 1g/4ml of water, yielding a clear, colorless solution
pKa:	6.1



AxyGen BioSciences
33210 Central Ave., Union City, CA 94587
Ph.: (510) 476-1090 Fax: (510) 476-1091
Email: sales@axygenbio.com
Website: www.axygenbio.com