



Certificate of Analysis

CP-Pre/Adi-65

AxyGen's ready to use buffers, precipitants, salts, additives, and detergents provide the most complete and economical solution to crystallization of proteins, peptides, small molecules and nucleic acids. They can be used as a rapid Initial screening method or to accelerate the creation of optimization/follow up screens. Reagents are supplied in sterile bottles (200ml/bottle).

Product Name:	5M NaCl
Description:	Sodium Chloride dissolved in water
Catalog #:	CP-Pre/Adi-65
Chemical Name:	Sodium Chloride
Molecular Formula:	NaCl
Molecular Weight:	58.44
CAS Number:	7647-14-5
Beilstein Registry:	3534976
EC Number:	2315983
Appearance:	White powder
Purity:	≥99.5%
MP:	801 °C
Solubility:	Dissolves at 1M of water, yielding a clear, colorless solution



AxyGen BioSciences
33210 Central Ave., Union City, CA 94587
Ph.: (510) 476-1090 Fax: (510) 476-1091
Email: sales@axygenbio.com
Website: www.axygenbio.com